

TECHNICAL DATA SHEET

1. TRADE NAME: **MOUSE STOP GLUE TRAP**

Banana flavoured monitoring device (possibility of hazelnut, cherry, cinnamon and chocolate flavor) for the following target pests: *Mus musculus musculus*, *Mus musculus domesticus* and other small rodents.

Monitoring device for the IPM application complying with UNI EN 16636:2015 Regulation, that must be applied to control regularly the pests presence in a specific place for a defined time.

2. PRODUCT CODE: **P-02001DIN**

3. PICTURE:



4. TECHNICAL SPECIFICATIONS:

Mouse Stop Glue Trap is a monitoring device (glue trap).

It is composed of:

- Cardboard;
- Banana flavored glue (possibility of hazelnut, cherry, cinnamon and chocolate flavor);
- White, black-printed silicon release paper, indicating the instructions for its detachment.
- Flow pack packs;
- Cardboard display.

Trap size (not assembled trap):

- **Length:** 254 mm
- **Width:** 120 mm

TECHNICAL DATA SHEET

Glue Quantity: 15 g

Glue area: 90 x 185 mm

Colour: As in the picture

5. HOW TO USE:

Place Mouse Stop Glue Trap in the following places:

- Near cable trough and electric boards, dishwashers, doors and slots, lofts, cellars, garages, etc.;
- Along the rooms perimeters to be disinfested and along the rodent habitual paths, which can be identified through excrements, fur, dens, paw prints, etc.
- Along the walls.

Theoretical covered area: at a distance of 5-10 meters the one from the other, according to the assumed pest infestation.

Suggested monitoring period: all the yearlong, especially in the period between October and April, in which it is registered peaks due to the external temperature reduction.

Suggested frequency of the controls: weekly. Change the trap after three months (at the latest) or before according to the pest infestation or the presence of dust, water and/or damp in the environment.

Storage conditions: in a cool and dry room, temperature below 25°C and R.H. below 55%.
If well stored (intact packaging), the product shelf life is 3 years.

Information on the Best Practice: Can be used in agri-food environments.

Evaluate the use of the product as it is or inside special containers, depending on the type of environment to be monitored.

The use inside containers allows to protect from dust and also to reduce the chances of catching non-target animals.

In order to respect animal welfare, we recommend a frequency of control of the adhesive traps equal to several times / day (at least 2), as recommended by the Guideline on Best Practice in the Use of Rodenticide Baits as Biocides in the European Union, issued in September 2013 from EBPF (European Biocidal Products Forum).

If, in the control phase of the adhesive trap, a mouse was found alive, it must be killed quickly and humanely. A strong blow to the head is the most advisable method according to the 2006 Animal Welfare Act UK. Absolutely avoid drowning.

If catches of non-target organisms are found, contact a wild animal recovery center (or a wildlife center). The removal of the live animal by unskilled personnel can cause further suffering and injury.



DATE: 04/12/2019 REV N°: 2 MOD: 4.1.3. UK

TECHNICAL DATA SHEET

Try to remove the animal with liquid paraffin, and then clean it with soap and water, only when it has been impossible to find professional help.

6. PACKAGING AND SPARE PARTS:

- 2 traps per flow pack;
- 24 packs per display.

7. DISPOSAL CONSIDERATIONS

Waste disposal should be in accordance with existing Community, National and local regulations.